


ORIHOMEDICAL

G m b H · I M P L A N T S



**DYNAMIC FIXATEUR EXTERNE
MODULAR SYSTEM**

Orthomedical axial movable head external fixators for femur and tibia, in adults



Cat. No.:	Type	Length mm min.	Length mm max.
96-410.10	short	230	250
96-410.20	medium	250	290
96-410.30	long	290	370

Orthomedical axial movable head external fixators for humerus, in pediatric



Cat. No.:	Type	Length mm min.	Length mm max.
96-411.10	short	190	210
96-411.20	medium	210	250
96-411.30	long	230	290

Orthomedical axial movable T-head external fixators for tibia, in adults



96-412.50
T-Clamp

Cat. No.:	Type	Length mm min.	Length mm max.
96-412.10	short	200	220
96-412.20	medium	220	260
96-412.30	long	260	340

Orthomedical axial movable T-head external fixators for humerus, in pediatric



96-413.50
T-Clamp



96-413.10
short



96-413.20
medium



96-413.30
long

Orthomedical axial movable head 135° external fixators for femur and tibia, in adults



Cat. No.:	Type	Length mm min.	Length mm max.
96-414.10	short	244	265
96-414.20	medium	264	305
96-414.30	long	304	385

Orthomedical axial movable square head external fixators for femur and tibia, in adults



Cat. No.:	Type	Length mm min.	Length mm max.
96-415.10	short	220	240
96-415.20	medium	240	280
96-415.30	long	280	360

Orthomedical limb lengthening external fixators for femur and tibia



Cat. No.:	Type	Length mm min.	Length mm max.
96-416.10	short	200	250
96-416.20	medium	240	330
96-416.30	long	280	410

Orthomedical limb reconstruction fixators



Cat. No.:	Type	Length mm min.	Length mm max.
96-417.10	short	300	320
96-417.20	medium	320	370
96-417.30	long	360	420

Orthomedical limb reconstruction pediatric fixators



96-418.10
short



96-418.20
medium



96-418.30
long

Orthomedical dynamic axial ankle fixators, in adults



96-419.10
left



96-419.20
right

Orthomedical external fixators for femur, valgus-varus




Cat. No.:	Type
96-420.10	short
96-420.20	long

Orthomedical external fixators for pelvic



96-421.10

Orthomedical external fixator complete set



96-450.00

Cat. No.:	Description		Pieces
96-410.10	Orthomedical dynamic axial moveable head external fixator for femur and tibia, in adults	short	1
96-410.20	Orthomedical dynamic axial moveable head external fixator for femur and tibia, in adults	medium	1
96-410.30	Orthomedical dynamic axial moveable head external fixator for femur and tibia, in adults	long	1
96-411.10	Orthomedical dynamic axial moveable head external fixator for humerus, in pediatric	short	1
96-411.20	Orthomedical dynamic axial moveable head external fixator for humerus, in pediatric	medium	1
96-411.30	Orthomedical dynamic axial moveable head external fixator for humerus, in pediatric	long	1
96-412.10	Orthomedical dynamic axial moveable T-head external fixator for tibia, in adults	short	1
96-414.30	Orthomedical dynamic axial moveable T-head 135° external fixator for femur and tibia, in adults	long	1
96-419.20	Orthomedical axial ankle fixator, in adults		
96-430.13	Schanz nail cortical screws for Orthomedical Fix	ø 3.5 x 30 x 130	6
96-430.15	Schanz nail cortical screws for Orthomedical Fix	ø 4.0 x 35 x 150	6
96-430.18	Schanz nail cortical screws for Orthomedical Fix	ø 4.0 x 40 x 180	6
96-430.20	Schanz nail cortical screws for Orthomedical Fix	ø 5.0 x 45 x 180	6
96-430.22	Schanz nail cortical screws for Orthomedical Fix	ø 6.0 x 45 x 180	6
96-431.18	Schanz nail cancellous screws for Orthomedical Fix	ø 6.0 x 45 x 180	6
96-433.01	Trochar	ø 6.0	1
96-433.05	Hexagonal Wrench	ø 5.0	1
96-433.06	Hexagonal Wrench	ø 6.0 long	1
96-433.02	Schanz T-Wrench		1
96-433.09	Drill bit with stop	ø 3.2	1
96-433.10	Drill bit with stop	ø 4.5	1
96-433.11	Drill	ø 2.5	1
96-433.03	T-Wrench (Quick Coupling Drill)		1
96-433.15	Drill Guide	ø 3.2 x 40	1
96-433.16	Drill Guide	ø 3.2 x 80	1
96-433.17	Drill Guide	ø 4.5 x 40	1
96-433.18	Drill Guide	ø 4.5 x 80	1
96-433.19	Sleeve Guide	ø 6.0 x 60	6
96-433.20	Sleeve Guide	ø 6.0 x 100	6

Schanz Screw with cortical thread



Cat. No.:	Length mm	Size ø
96-430.13	30 x 130	3.5
96-430.15	35 x 150	4.0
96-430.18	40 x 180	4.0
96-430.20	45 x 180	5.0
96-430.22	45 x 180	6.0

Schanz Screw with cancellous thread



Cat. No.:	Length mm	Size ø
96-431.18	45 x 180	6.0



96-433.01 Trochar



96-433.02 Schanz T-Wrench



96-433.05 Hexagonal Wrench 5.0 mm
96-433.06 Hexagonal Wrench 6.0 mm



Drill Bit with stop

Cat. No.:	Size ø
96-433.09	3.2
96-433.10	4.5



Drill

Cat. No.:	Size ø	Length mm
96-433.11	2.5	110 x 50



Drill Guide

Cat. No.:	Size ø	Length mm
96-433.15	3.2	40
96-433.16	3.2	80
96-433.17	4.5	40
96-433.18	4.5	80



Sleeve Guide

Cat. No.:	Size ø	Length mm
96-433.19	6.0	60
96-433.20	6.0	100



96-433.03 T-Wrech (Quick Coupling Drill)

ORIHOMEDICAL

G m b H · I M P L A N T S

Hauptstrasse 5 · D-78589 Dürbheim
Tel.: 0 74 24 / 94 03-40 · Fax: 0 74 24 / 94 03-43
www.orthomedical.de · email: mail@orthomedical.de
